



### Protocol for Removing PPE:

Wash hands and remove visible bodily fluids/blood contaminants on the outer surfaces of both hands → Wash hands replace outer gloves with new gloves → Remove powered air-purifying respirator or self-priming filter-type full-face mask/mask (if used) → Wash hands → Remove disposable gowns along with outer gloves (if used) → Wash hands and put on outer gloves → Enter Removal Area No. ① → Wash hands and remove protective clothing along with outer gloves (for gloves and protective clothing, turn inside out, while rolling them down) (note: if used, remove the waterproof boot covers with clothing) → Wash hands → Enter Removal Area No. ② → Wash hands and remove goggles → Wash hands and remove mask → Wash hands and remove cap → Wash hands and remove inner disposable latex gloves → Wash hands and leave Removal Area No. ② → Wash hands, take a shower, put on clean clothes and enter the clean area